



DEMOGRAPHICS

PRISCILLA
FEMALE
43 YEARS OLD

DIAGNOSIS

CLASS III
ANTERIOR CROSSBITE
BILATERAL POSTERIOR
CROSSBITE

MECHANICS

INVISALIGN®
51 TOTAL ALIGNERS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
24-36 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS

DEMOGRAPHICS

PRISCILLA
FEMALE
43 YEARS OLD

DIAGNOSIS

CLASS III
ANTERIOR CROSSBITE
BILATERAL POSTERIOR
CROSSBITE

MECHANICS

INVISALIGN[®]
51 TOTAL ALIGNERS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
24-36 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDDENT:
12 MONTHS



BASELINE



DEMOGRAPHICS

PRISCILLA
FEMALE
43 YEARS OLD

DIAGNOSIS

CLASS III
ANTERIOR CROSSBITE
BILATERAL POSTERIOR
CROSSBITE

MECHANICS

INVISALIGN[®]
51 TOTAL ALIGNERS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
24-36 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDDENT:
12 MONTHS

PROGRESS

6 weeks with first
6 aligners



DEMOGRAPHICS

PRISCILLA
FEMALE
43 YEARS OLD

DIAGNOSIS

CLASS III
ANTERIOR CROSSBITE
BILATERAL POSTERIOR
CROSSBITE

MECHANICS

INVISALIGN[®]
51 TOTAL ALIGNERS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
24-36 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS



PROGRESS Pre-Refinement



DEMOGRAPHICS

PRISCILLA
FEMALE
43 YEARS OLD

DIAGNOSIS

CLASS III
ANTERIOR CROSSBITE
BILATERAL POSTERIOR
CROSSBITE

MECHANICS

INVISALIGN[®]
51 TOTAL ALIGNERS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
24-36 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS



FINISH

