

CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON® SELF-
LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT TIME:
21 MONTHS

ACTUAL TREATMENT TIME WITH
ACCELEDENT:
8 MONTHS

“He was very compliant in all aspects of treatment. He was very excited about AcceleDent and due to its use was able to progress very efficiently through his orthodontic treatment.”

- Dr. Jason Kaplan

CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
8 MONTHS



BASELINE

Chief Concern: Cosmetic needs and making his smile more aesthetic



CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
8 MONTHS



FINISH

Remarks: *He was very compliant in all aspects of treatment. He was very excited about AcceleDent and due to its use was able to progress very efficiently through his orthodontic treatment.*



CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
8 MONTHS



RESULTS: 8 MONTHS

CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
8 MONTHS



RESULTS: 8 MONTHS

CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
8 MONTHS



RESULTS: 8 MONTHS

CASE STUDY



DEMOGRAPHICS

MALE
32 YEARS OLD

DIAGNOSIS

CLASS I
CROWDING
ANTERIOR CROSSBITE
POSTERIOR CROSSBITE

MECHANICS

INSIGNIA™ DAMON®
SELF-LIGATING BRACKETS

TREATMENT TIME

PROJECTED TREATMENT
TIME:
21 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
8 MONTHS

